



REDISCOVER COMFORT IN LIFE



ABOUT US

Dynamic Techno Medicals is one of the largest manufacturers & marketers of medical devices in India with a well diversified product portfolio that includes orthopaedic appliances, surgical dressings, sports supports, rehabilitation aids, medical compression stockings, silicone foot care products, medical disposables and unique consumer care products

OTHER SALIENT FEATURES:

- Started operations in 1980
- Headquartered in Kerala, India, close to Cochin International Airport
- 10 state-of-the-art manufacturing facilities
- 1,15,000 square feet of production floor space
- 20,000 square feet of warehousing space at factories
- Class 10,000 and 100,000 clean rooms
- Exporting to more than 35 countries, spanning 6 continents
- Supplying under private label for major corporations in Europe and Australia

ISO 9001: 2008, 13485: 2003 Company | Most products are  certified

Large R & D dept, company spends 5% of its total turnover on R & D activities
Dynamic has been started on a promise. An assurance that we will do whatever it takes to make quality health care affordable and accessible to all. In pursuit of this, we create unique products to help people

“ REDISCOVER COMFORT IN LIFE”.

GLOBAL PRESENCE



DYNAMIC)
TECHNO MEDICALS

Kodikuthimala, Aluva, Kerala
Phone: +91 484 2837970, 2837788
info@dynamictechnomedicals.com



CONTENTS

04 3D KNITTED PRODUCTS

08 NECK

12 CHEST & SHOULDER

18 ABDOMEN

20 BACK

24 ARM & WRIST

34 KNEE & ANKLE

40 SURGICAL DRESSINGS & BANDAGES

44 ADVANCED DRESSINGS

50 COMPRESSION STOCKINGS

56 TRACTION APPLIANCE

60 REHABILITATION AIDS

64 PODIATRY

68 MOTHER CARE

74 SPORTS SUPPORTS

82 MISCELLANEOUS



3D KNITTED PRODUCTS



1225

Sizes: Small, Medium, Large, X-Large, XX-Large
XXX-Large

GENUGRIP® KNEE BRACE

- 3D knitted anatomically contoured breathable brace for adequate compression
- Anatomically shaped silicone pressure pad for pain relief
- Soft knitting at the popliteal region for excellent wearing comfort
- Stitch-free, fused pocket makes the brace appealing
- Ideal for sprain, strain, patellar dislocation, arthritis and post operation

Left & Right



1222

Sizes: Small, Medium, Large, X-Large, XX-Large
XXX-Large

GENUGRIP® ML KNEE BRACE WITH MEDIO LATERAL STAYS

- 3D knitted, anatomically contoured brace for effective compression
- Anatomically shaped silicone pressure pad for pain relief
- Conforms to body contour for secure fit and wearing comfort
- Medio lateral spiral stays provide MCL/ACL support



1228

Sizes: Small, Medium, Large, X-Large, XX-Large
XXX-Large

GENUGRIP® HB HINGED KNEE BRACE

- 3D knitted, four way stretchable material prevents constriction, offers adequate compression and increases blood circulation
- Breathable knit provides air circulation and reduces sweat accumulation
- Medio lateral metal hinges mimics the natural knee joint to support the cruciate ligament during flexion
- Opposing straps in the anterior and posterior region provides reinforcement and joint stability

Left & Right



1323

Sizes: Small, Medium, Large, X-Large, XX-Large
XXX-Large, XXXX-Large

LUMBOGRIP® LUMBO SACRAL CORSET

- 3D knitted, breathable, stretchable brace for adequate compression
- Posterior metal stays facilitate correct posture
- Fastening handles enable easy application
- Ideal for chronic low back pain, Osteoporotic pain, Post operative rehabilitation and discectomy



1363

Sizes: Small, Medium, Large, X-Large, XX-Large
XXX-Large, XXXX-Large

LUMBOGRIP® SL LUMBAR BRACE WITH SILICONE PAD

- 3D knitted anatomically contoured breathable brace for adequate compression
- Anatomically shaped silicone pressure pad for pain relief
- Plastic stays on either sides provide additional support and rigidity
- Ideal for chronic low back pain, Osteoporotic pain, Post operative rehabilitation and discectomy



1365

Sizes: Small, Medium, Large, X-Large, XX-Large
XXX-Large, XXXX- Large

LUMBOGRIP® DS LUMBO SACRAL CORSET (DUAL STRAPS)

- 3D knitted, anatomically contoured brace for effective compression
- The double reinforced straps allow individually adjustable compression to support core muscles and to ensure an adaptive fit for better compliance
- Anatomically contoured posterior metal stays facilitate correct posture and permit individual adaptation to the patient's shape
- Breathable elastic provides air circulation and reduces sweat accumulation



1676

Sizes: Small, Medium, Large, X-Large, XX-Large

EPIGRIP® ELBOW BRACE

- 3D knitted anatomically contoured breathable brace for adequate compression
- Anatomically shaped silicone pressure pad for pain relief
- Stitch-free, fused pocket makes the brace appealing
- Ideal for osteoarthritis, sprain, strain, post traumatic and operative corrections



1207

Sizes: Small, Medium, Large, X-Large, XX-Large

MALLEOGRIP® ANKLE BRACE

- 3D knitted anatomically contoured breathable brace for adequate compression
- Two anatomically shaped silicone pressure pads for pain relief
- Stitch-free, fused pocket makes the brace appealing
- Ideal for ligament, soft tissue and tendon injuries and as prophylaxis during sports activity



1209

Sizes: Small, Medium, Large, X-Large, XX-Large

MALLEOGRIP'S® ANKLE BINDER

- 3D knitted, anatomically contoured brace for effective compression
- Breathable elastic strap wraps in a "figure of 8" pattern
- The strap system stabilizes the ankle

NECK



1001

Size: Paediatric, Child, Small, Medium, Large, X-Large

DYNĀ SOFT CERVICAL COLLAR

- Designed for minimal immobilisation during sprain or strain
- Made of high-density Polyurethane foam with Polyethylene sheet reinforcement
- Plastic bushes for air circulation



1003

Size: Small, Medium, Large, X-Large

DYNĀ INNO-LIFE SOFT CERVICAL COLLAR

- High density moulded EVA foam padded and covered with cotton stockinette
- Removable specially knitted tubular cover can be washed, thereby promoting patient hygiene
- Designed for minimal immobilisation during sprain or strain



1002

Size: Small, Medium, Large, X-Large



SILVER SOFT CERVICAL COLLAR

Silver collar designed for minimal immobilisation

THE SILVER EFFECT

- Anti odour** : Promotes patient hygiene by eliminating odour
- Anti microbial** : Eliminates 99.9% bacteria in less than 1 hour
- Anti itch** : Protects the skin against infection and rashes
- All natural** : Silver is safe and non-toxic, contain no chemicals
- Permanent** : The effect of silver does not decrease on washing or wearing off





1010

Size: Small, Medium, Large, X-Large

DYNÁ[®] CONTOUR CERVICAL COLLAR

- Maintains head in neutral position during travelling or prolonged sitting
- Made of firm foam with cotton stockinette covering



1050

Size: Paediatric, Child, Small, Medium, Large, X-Large

DYNÁ[®] HARD CERVICAL COLLAR

- Made of high density Polyethylene sheet for moderate immobilisation
- Comes in two-piece design for height adjustment
- Ideal for moderate to severe cervical sprain



1090

Size: Small, Medium, Large, X-Large

TOPPHIL[®] CERVICAL IMMOBILISER

- Two-Piece collar designed for optimal immobilisation
- X-ray lucent and anterior tracheal opening for easy access to carotid pulse

TOPPHIL[®] Restricts

Flexion/Extension	Rotation	Lateral Flexion
71.1%	56.3%	33.6%



1150

Size: One size fits most

DYNÁ[®] AMBULANCE CERVICAL IMMOBILISER

- Single piece design aids easy wearing of the collar without disturbing the cervical region
- Ideal for prompt and temporary cervical immobilisation of accident victim during transportation



1080

Size: Small, Medium, Large, X-Large

DYNÁ® TRAUMA CERVICAL IMMOBILISER

- Ideal for post trauma cervical immobilisation
- Provides motion restriction post traumatic injury and surgical care
- Oversized tracheal opening for easy assessment



1004

Size: Child, Small, Medium, Large, X-Large

DYNÁ® CERVICAL COLLAR WRAP ROUND

- Anatomically contoured collar designed for semi-rigid support
- Supports the neck and prevents acute pain during travelling or prolonged sitting (as in front of a computer)



1100

Size: Universal

RELAXZON® THERAPEUTIC SUPPORT PILLOW

- Made of hypo-allergenic, high quality polyurethane foam with 100% cotton cover
- Provides proper sleeping posture by conforming to the natural sleeping curve of the neck
- Provides relief from stress related headaches




1020

Size: Universal

DYNÁ® AIRWAY COLLAR

- Single piece design aids easy wearing of the collar without disturbing the cervical region
- Ideal for prompt and temporary cervical immobilisation of accident victim during transportation

CHEST AND SHOULDER




1401

Size: Child, Paediatric, Small, Medium, Large, X-Large, XX-Large

DYNA[®] CLAVICLE BRACE

- Bilateral underarm axillary pad cushioning to prevent excess pressure
- Ideal for the correction and management of clavicular fractures, injuries & stooped shoulders



1405

Size: Child, Small, Medium, Large, X-Large, XX-Large

DYNA[®] INNO-LIFE CLAVICLE BRACE

- Bilateral underarm axillary pad cushioning to prevent excess pressure
- Fully padded foam straps for maximum patient comfort
- Ideal for the correction and management of clavicular fractures, injuries & stooped shoulders



2800

Size: Small, Medium, Large, X-Large, XX-Large

SEGO[®] RIB BRACE PLAIN

- Made of 8" width high quality, soft and non-irritating elastic material
- Offers adequate compression and support to the ribs
- Ideal for post fracture stabilisation of the sternum

Available in male & female models



2820(M)
2830(F)

Size: Small, Medium, Large, X-large, XX-Large

SEGO® BREATH RIB BRACE PLAIN

- Made of breathable elastic material with air vents allowing air circulation
- Offers adequate compression and support to the ribs
- Ideal for post fracture stabilisation of the sternum

Available in male and female models



1435

Size: Small, Medium, Large, X-large, XX-Large, XXX-Large

DYNA® CHEST BRACE WITH STERNAL PAD

- Full elastic brace with sternal pad made of soft foam to hold the rib brace intact
- Prevents strain on the chest while coughing, sneezing or deep breathing




1436

Size: Small, Medium, Large, X-Large, XX-Large, XXX-Large

DYNA® CHEST BRACE WITH SILVER STERNAL PAD

- Full elastic brace with sternal pad made of soft foam to hold the rib brace intact
- Sternal pad is covered with anti-microbial silver fabric
- Silver reduces the chance of infection on the operated wound
- Prevents strain on the chest while coughing, sneezing or deep breathing



1410

Size: Child, Small, Medium, Large, X-Large

DYNÁ[®] THORACO LUMBAR APPLIANCE

- Made of light weight Polypropylene for better stability of the spine
- The brace extends from thoracic to lumbosacral region for optimal spinal support
- Anatomically contoured and well padded for maximum patient comfort
- Maintains spine in neutral position and permits ease of movement



1411

Size: Child, Small, Medium, Large, X-Large

DYNÁ[®] THORACO LUMBAR APPLIANCE NOVA

- Made of manmade fabric reinforced with light weight & durable Aluminium
- The brace extends from thoracic to lumbosacral region for optimal spinal support
- Anatomically contoured and well padded for maximum patient comfort
- Maintains spine in neutral position and permits ease of movement



1412

Size: Small, Medium, Large, X-large, XX-Large

DYNÁ[®] BREATH THORACO LUMBAR APPLIANCE

- Made of breathable elastic material reinforced with light weight & durable Aluminium
- The brace extends from thoracic to lumbosacral region for optimal spinal support
- 2 posterior and 2 lateral Aluminium bars maintain spine in neutral position and permits ease of movement



1420

Size: Child, Small, Medium, Large, X-Large

DYNA[®] SPINAL JACKET

- Made of light weight polyethylene for better stability of a spine
- Anatomically contoured and well padded for maximum patient comfort
- The anterior and posterior support provides maximum immobilisation



1441

Size: One size fits most

DYNA[®] HYPER EXTENSION BRACE NEW

- Two lateral pads for trunk stability
- Hinged sternal and pubic pads adjustable both horizontally and vertically
- Made of aluminium, light weight and comfortable to use
- Ideal for any problem which needs moderate immobilisation of lower thoracic to upper lumbar region



1445

Size: Child, Small, Medium, Large, X-Large

DYNA[®] INNO- LIFE HYPER EXTENSION BRACE

- Made of light weight Aluminium, with soft sternal, pubic and lateral pads
- Wide and soft spinal pad for better comfort
- The product can be adjusted both vertically & horizontally as per patient's requirement
- Ideal for stable compression fractures and post operative immobilisation of the spine and adolescent kyphosis

1443

Size: One size fits most



DYNÁ INNO- LIFE HYPER EXTENSION BRACE PREMIUM

- Adjustable, light weight and comfortable brace for optimal spinal support
- Soft sternal pubic & lateral pads with wide & soft posterior spinal pad
- Single side closing mechanism for user friendly application



ABDOMEN



2700

Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

Cling® POST MATERNITY CORSET 25

- Made of premium quality elastic material
- 25cm wide corset covers the entire abdomen
- Supports weakened abdominal muscles post surgery/delivery



2720

Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

Cling® BREATH POST MATERNITY CORSET

- Made of premium quality breathable elastic material
- Special wide weave aids air circulation and prevents body odour
- 25cm wide corset covers the entire abdomen
- Supports weakened abdominal muscles post surgery/delivery



2730

Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

SEGO® ABDOMINAL CORSET PLAIN

- 20cm wide corset made of premium elastic material
- Fits the body contours and tightens abdominal muscles



2710

Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

SEGO® BREATH ABDOMINAL CORSET PLAIN

- 20cm wide corset made of breathable elastic material
- Fits the body contours and tightens abdominal muscles



2734

Sizes: Universal

SEGO® ELASTIC ABDOMINAL CORSET (4 PANEL)

- Hook compatible elastic permits a closure at any point on the circumference of the binder
- Four panel design help in vertical size adjustments and reduce rolling
- Made of full elastic, soft material to provide comfort and flexibility
- Designed for use by both men and women



BACK



1320


Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

Elnova® SURGICAL LUMBO SACRAL CORSET

- Full elastic corset fits snugly around the body contours
- Four anatomically contoured metal stays can be individually adapted to suit the natural curvature of the patient's spine
- Back pad with increased length covers the entire LS region for better support and stability
- Additional strap enhances proper fit
- Ideal for chronic low back pain and lumbar disc problems



1321



1320

Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large, XXXX-Large

Elnova® SURGICAL LUMBO SACRAL CORSET BROADER

- Full elastic corset fits snugly around the body contours
- Four anatomically contoured metal stays can be individually adapted to suit the natural curvature of the patient's spine
- Back pad with increased width covers the entire LS region for better support and stability
- Additional strap enhances proper fit
- Ideal for chronic low back pain and lumbar disc problems





1301

Sizes: 28" to 58" in 2" intervals

DYNA[®] SURGICAL LUMBO SACRAL CORSET

- Made of sanforised cotton with elastic in the front
- Two rigid and semi-rigid posterior bars for spinal stability



1331

Sizes: Small, Medium, Large, X-Large, XX-Large
XXX-Large

DYNA[®] BREATH X SURGICAL LUMBO SACRAL CORSET

- Made of premium quality breathable elastic material
- Special wide weave aids air circulation and prevents body odour
- Criss cross straps at the LS region and two rigid posterior bars at either sides of the spine provides additional support
- For chronic low back pain, strain, osteoporotic pain, disc diseases etc.



1340

Sizes: Small, Medium, Large, X-Large, XX-Large

DYNA[®] INDUSTRIAL LS BELT

- Made of high quality elastic with contour fit at the lumbo sacral region
- 2 semi rigid posterior bars for rigidity and to prevent exaggerated movements
- 2 shoulder straps for suspension and maintaining the position
- Ideal for labourers in jobs involving frequent or heavy lifting, carrying, pulling and bending



1301

Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

DYNA[®] SURGICAL LUMBO SACRAL CORSET

- Elastic supports both sides for snug fit
- Front and back panels are made of soft foam for comfort
- Hook & loop closures for easy application and removal
- Two rigid aluminium bars and two plastic bars for reinforcement

1328

Sizes: One size fits most

DYNA® PAVLIK HARNESS(CDH SPLINT)

- Soft padded straps avoid irritation and gives comfort to the baby
- For Congenital Dysplasia of the Hip of children up to 8 months



1338

Sizes: One size fits most

DYNA® UMBILICAL TRUSS

- Reduces hernia in children by giving an external support to the weakened muscles
- Soft material keeps the baby's body warm
- Hook & loop straps for easy application



1335

Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

DYNA® SURGICAL TRUSS

- Flexible compression pads with tubular knit perineal straps
- Gentle aid for inguinal hernia



1342

Sizes: Universal

DYNA® ORTHOPAEDIC HEATING BELT

- With cover and belt for hands - free operation
- Ultra soft washable cotton cover for user comfort
- Useful for backaches, sprain, muscular and joint pains, abscesses, Rheumatoid or other arthritis, geriatric care, sore over worked muscles





ARM AND WRIST



1605

Size: Paediatric, Child, Adult, Small, Medium, Large, X-Large

DYNÁ[®] ARM SLING DELUXE

- Made of special cloth material with adjustable shoulder straps
- Thumb loop holds the forearm in position
- Ideal support for the forearm and wrist following soft tissue, bone and ligament injuries and fractures
- Effective aid for conservative treatment of forearm fractures



1601

Size: Child, Small, Medium, Large, X-Large

DYNÁ[®] ARM SLING

- Traditional sling made of sanforised cotton
- Additional support for patients with forearm plaster-of-paris cast
- Ideal support for the forearm and wrist following soft tissue, bone and ligament injuries and fractures



1601

Size: Child, Small, Medium, Large, X-Large

DYNÁ[®] ARM SLING

- Traditional envelope arm sling made in soft foam
- Additional support for patients with forearm plaster-of-paris cast
- Ideal support for the forearm and wrist following soft tissue, bone and ligament injuries and fractures



1610

Size: Child, Small, Medium, Large, X-Large

DYNÁ[®] ARM SLING POUCH

- Contoured hand rest area and thumb loop hold the forearm in position
- Well-padded adjustable shoulder strap reduces stress and soreness to the neck
- Ideal support for the forearm and wrist following soft tissue, bone and ligament injuries and fractures



1635

Size: One size fits most

DYNA® SLING WITH TIE

- Soft foam in sleeved stockinette with plastic fastening ties
- Used as a traditional sling to provide effective support without the need for additional padding



1637

Size: 5cm Width X 6m Length X 2

DYNA® CUFF - N - COLLAR

- Light- weight polyurethane foam in a stockinette for comfortable support
- Available in 6m rolls in a dispensing carton with 5 plastic fastening ties



1625

Size: Small, Medium, Large, X-Large, XX-Large

DYNA® SHOULDER IMMOBILISER

- Traditional sling designed to hold the arm in 90° angle
- Made of sanforised cotton with shoulder straps
- Circumferential arm wrap for complete immobilisation of the shoulder



1626

Size: Small, Medium, Large, X-Large, XX-Large

DYNA® SHOULDER IMMOBILISER SPECIAL

- An ideal aid to immobilise the shoulder, post shoulder dislocation, subluxation and post operative conditions
- 2 elastic straps attached to the thoracic brace helps to hold the hand in 90° angle
- Special design removes strain on the neck



1623

Size: Small, Medium, Large, X-Large, XX-Large

DYNA[®] ELASTIC SHOULDER IMMOBILISER

- Elastic for comfort and snug fit
- Arm straps keep hand in position
- Additional shoulder cup ensures complete immobilisation



1624

Size: Universal

DYNA[®] SHOULDER IMMOBILISER

- For easy application and removal
- Made of soft foam for comfort
- Long soft wrap for maximum immobilisation
- Universal design with hook & loop adjustments



1628

Size: Small, Medium, Large, X-Large, XX-Large

DYNA[®] BREATH SHOULDER IMMOBILISER SPECIAL

- Breathable elastic brace provides air circulation and prevents body odour
- An ideal aid to immobilise the shoulder, post shoulder dislocation, subluxation and post operative conditions
- 2 elastic straps attached to the thoracic brace helps to hold the hand in 90° angle
- Special design removes strain on the neck



2950

Size: Small, Medium, Large, X-Large

SEGO[®] SHOULDER SUPPORT

- Anatomically shaped 4-way stretchable elastic fabric
- Wide elasticated band ensures easy application
- Provides support, compression and therapeutic warmth to the shoulder and aids in healing of the shoulder injuries
- Minimises the chances of further shoulder injuries



2900

Size: Small, Medium, Large, X-Large, XX-Large

SEGO® ELBOW SUPPORT

- 4-way stretchable brace prevents constriction and improves blood circulation
- Unique single way knit at the elbow joint ensures free and easy movement
- Anatomically contoured design prevents slipping and bunching
- Ideal for sprain, strain, arthritis and weak and over stressed elbow



1670

Size: One size fits most

DYNA® TENNIS ELBOW BRACE

- Fits around the forearm and takes away strain from the lateral epicondyle
- Suitable for all sporting activities for both men and women
- Can be adjusted to fit most male and female athletes
- Comfortable and durable



1671

Size: One size fits most

DYNA® TENNIS ELBOW BRACE X-TRA

- Fits around the forearm and takes away strain from the lateral epicondyle
- Extra width offers extra support to the forearm
- Suitable for all sporting activities for both men and women
- Can be adjusted to fit most male and female athletes



1672

Size: One size fits most

DYNA® INNO-LIFE TENNIS ELBOW BRACE X-TRA

- Made of breathable material that aids air circulation and prevents body odour
- Fits around the forearm and takes away strain from the lateral epicondyle
- Extra width offers extra support to the forearm
- Suitable for all sporting activities for both men and women



1673

Size: 1, 2

DYNA[®] INNO-LIFE TENNIS ELBOW BRACE

- Provides adequate compression to help relieve tennis elbow pain and discomfort
- Interior shock absorbing cushion targets compression exactly where needed to help provide relief
- Suitable for all sporting activities for both men and women
- Adjustable fasteners for customised fit



2920

Size: Small, Medium, Large, X-Large, XX-Large

SEGO[®] WRIST SUPPORT

- Unique "Pull and Wrap" design ensures proper fit
- Ideal for tendon injuries and arthritic conditions of the wrist



2918

Size: One size fits most

SEGO[®] WRIST WRAP WITH THUMB LOOP

- Rigid wrap around design provides customised compression and fitting
- Anatomical thumb opening provides slight abduction for better comfort of the thumb
- Ideal for strain, sprain and chronically weak wrist

Available in Black & Beige Colour



1640

Size: Small, Medium, Large, X-Large

DYNA[®] WRIST BRACE

- Made of superior elastic with anatomically contoured palmar stays reinforcement
- Maintains wrist in its functional position
- For Wrist Drop, Carpal Tunnel Syndrome, Colles' fracture and Arthritis

Left & Right



1642

Size: I,II,III

DYNA® WRIST BRACE REVERSIBLE

- Made of full elastic, with anatomically contoured rigid palmar stays
- Elastic wrist wrap provides added reinforcement
- Can be worn on both hands
- For Wrist Drop, Carpal Tunnel Syndrome, Colles' fracture and arthritis



1644

Size: Small, Medium, Large, X-Large

DYNA® BREATH WRIST BRACE

- Breathable elastic brace provides air circulation
- Made of superior elastic with anatomically contoured palmar stays reinforcement
- Maintains wrist in its functional position
- For Wrist Drop, Carpal Tunnel Syndrome, Colles' fracture and Arthritis

Left & Right



1645

Size: Small, Medium, Large, X-Large

DYNA® INNO-LIFE WRIST BRACE

- Made of superior elastic with anatomically contoured palmar stays reinforcement
- The palmar stay can be moulded to adapt to the contours of the patient's palm
- Maintains wrist in its functional position
- Ideal for Wrist Drop, Carpal Tunnel Syndrome, Colles' fracture and Arthritis

Left & Right



1654

Size: X- Small, Small, Medium, Large, X-Large

DYNA® INNO-LIFE BREATHABLE WRIST BRACE

- Breathable fabric aids air circulation and prevents body odour
- Removable metal palmar and dorsal stays for better support and complete immobilisation
- Additional wrist wrap for added reinforcement
- Palmar and dorsal stays allow minor customisation and accommodate larger arms

Left & Right



1660

Size: Small, Medium, Large, X-Large

DYNA[®] COCK UP BRACE

- Plastic coated malleable Aluminium with foam padding keeps the wrist in extension and provides comfort
- To correct hand and finger fractures, injuries and other deformities



1661

Size: Small, Medium, Large, X-Large

DYNA[®] DYNAMIC COCK UP BRACE

- Made of well padded, unbreakable, washable and reusable material with optional finger-thumb attachments
- Spring system helps to adjust wrist tension which is ideal for physiotherapy after extensor tendon repair
- Effective for post fracture rehabilitation, burns to prevent contractures and for cerebral palsy patients



1662

Size: Small, Medium, Large, X-Large

DYNA[®] DYNAMIC COCK UP BRACE (WITH FINGER EXTENSION ASSIST)

- Made of well padded, unbreakable, washable and reusable material with optional finger-thumb attachments
- The detachable finger extension assist has spring system to adjust wrist tension
- Effective for physiotherapy, post fracture rehabilitation, burns to prevent contractures and for cerebral palsy patients

Left & Right



1688

Size: Small, Medium, Large

DYNA[®] KNUCKLE BENDER BRACE

- Detachable spring to adjust tension as per patient requirement
- Ideal for post fracture rehabilitation, to restore full range of motion and eliminate joint stiffness



1655

Size: One size fits most

DYNA[®] THUMB SPICA BRACE

- Immobilises the thumb and wrist while allowing other fingers to move
- Provide support for thumb injuries, reduce pain and inflammation and prevents further injuries



1657

Size: One size fits most

DYNA[®] INNO-LIFE THUMB SPICA BRACE

- Immobilises the thumb without inhibiting hand movement
- Provide support for thumb injuries, reduce pain and inflammation and prevents further injuries



1680

Size: One size fits most

DYNA[®] FINGER EXTENSION BRACE

- Plastic coated malleable Aluminium with polyfoam padding
- For immobilization of finger following fractures and injury of Proximal Intra Phalanges



1682

Size: One size fits most

DYNA[®] FINGER COT

- Plastic coated malleable Aluminium with polyfoam padding
- For collateral ligament injury and fracture of Distal Intra Phalanges



1683

Size: Small, Medium, Large

DYNA® FROG SPLINT

- Plastic coated malleable Aluminium with protective foam padding
- Unique 4 legged construction allows self application
- For mallet finger, swan neck deformity and problems in DIP joint



1685

Size: 2, 3, 4, 5.5, 6, 7

DYNA® STACK FINGER SPLINT

- Made of sturdy plastic for support and immobilisation
- Light-weight and perforated for patient comfort
- Ideal for immobilisation of finger following fracture and injuries



1687

Size: Small, Medium, Large

DYNA® FINGER PROTECTOR SPLINT

- Made of malleable Aluminium with foam padding for comfort
- Provides complete immobilisation and protection for the finger
- Ideal for treatment of finger fractures, sprains, dislocations or as protection following injury



1632

Size: Child, Small, Medium, Large, X-Large

DYNA® ELBOW IMMOBILISER

- Made of soft foam in a wrap-around design for optimal immobilisation
- Medio lateral and posterior bars for stability of the elbow joint
- For complete elbow immobilisation post injury and surgery

KNEE & ANKLE



2565

Size: X-Small, Small, Medium, Large, X-Large, XX-Large, XXX-Large

SEGO® KNEE SUPPORT PLAIN

- 4-way stretchable brace prevents constriction and improves blood circulation
- Unique single way knit at the knee joint ensures free and easy movement
- Anatomically contoured design prevents slipping and bunching
- Ideal for sprain, strain, joint effusions and pain due to arthritis



2563

Size: Small, Medium, Large, X-Large, XX-Large

SEGO® KNEE SUPPORT (OP)

- 4-way stretchable brace prevents constriction and improves blood circulation
- Anterior patellar opening relieves pressure on patella
- Soft circular cushion around the open patella gives extra padded protection to knee joint
- Ideal for sprain, strain, joint effusions and pain due to arthritis



1270

Size: X-Small, Small, Medium, Large, X-Large, XX-Large, XXX-Large

DYNÁ® HINGED KNEE BRACE (WITH PATELLA SUPPORT)

GENU ORTHO

- Medio lateral metal hinges and soft buttress on the patellar region
- Ideal for patellar dislocation or subluxation, tendonitis, OSD and chondro malacia patella



1271

Size: X-Small, Small, Medium, Large, X-Large, XX-Large, XXX-Large

DYNÁ® KNEE BRACE (MEDIO LATERAL STABILISER)

GENU ML

- Anatomically contoured brace with medio lateral spiral stays
- Ideal for medio lateral ligament support, meniscus damage, joint effusion



1232

Size: One size fits most

DYNA[®] LIMITED MOTION KNEE BRACE PREMIUM

- Allows Range of Motion (ROM) at 10° increments
- Provides desired ROM with individual angle adjustments
- For controlled immobilisation after ligament injuries, knee dislocations, joint Arthroplasty etc



1233

Size: One size fits most

DYNA[®] LIMITED MOTION KNEE BRACE SHORT

- Universal knee brace with flexion/extension lock
- Allows Range of Motion and complete immobilisation
- Adjustable thigh & calf bands and medio-lateral support for complete immobilisation
- Ideal for post-operative immobilisation, ACL/PCL/MCL & LCL injuries, knee dislocation and joint arthroplasty



1250

Size: One size fits most

DYNA[®] KNEE IMMOBILISER DELUXE

- Unique 3 panel universal style
- Made of soft fabric
- For complete immobilisation of the knee after injury or surgery

Product length :35cm, 50cm, 60cm



1240

Size: Child, Small, Medium, Large, X-Large, XX-Large

DYNA[®] KNEE IMMOBILISER

- Made of soft foam with medio lateral and anatomically contoured posterior bars
- Provides complete immobilisation of the knee in full extension



1242

Size: Child, Small, Medium, Large, X-Large, XX-Large
XXX-Large, XXXX-Large

DYNA[®] INNO-LIFE KNEE IMMOBILISER LONG

- Anatomically contoured posterior bars and medio-lateral bars offer complete immobilisation
- Patellar opening releases patellar pressure ensuring patient comfort
- Elastic patellar strap holds the patella in position and ensures immobilisation of the knee joint
- Extra length for complete immobilisation and added support

Available in Beige & Grey



1255

Size: Small, Medium, Large, X-Large, XX-Large
XXX-Large

DYNA[®] HINGED KNEE BRACE OPEN PATELLA

- Open patella design with soft circular cushion allows flexibility around the knee cap
- Medio lateral hinges ensure proper support and free movement
- Ideal for mild to moderate medio lateral instability of the knee, patellar dislocation, subluxation etc.

Available in Beige & Grey



1260

Size: Small, Medium, Large, X-Large, XX-Large
XXX-Large

DYNA[®] INNO-LIFE HINGED KNEE BRACE OPEN PATELLA

- 4 - way stretchable and breathable fabric provides comfort and compression
- Medio lateral hinges ensure proper support and free movement
- Open patella design for patellar tracking
- Wrap around design with anterior and posterior strapping ensures comfortable fit



1202

Size: Small, Medium, Large, X-Large

DYNA[®] INNO-LIFE ANKLE IMMOBILISER WITH LACE

- Anatomically designed, rigid splint ensures immobilisation of the ankle joint
- Lace design for easy and comfortable fit
- Sleek, light weight and bilateral design, can be used inside the shoe for both right & left leg



1201

Size: One size fits most

DYNA[®] ANKLE IMMOBILISER

- Solid outer shell provides complete immobilisation to the ankle
- Inside foam padding for comfort and better stability
- Ideal for inversion & eversion of foot and soft tissue injuries of the ankle



2509

Size: Small, Medium, Large, X-Large

SEGO[®] ANKLE BINDER

- The sleeve is made of sturdy, elastic fabric which provides warmth, mild compression and improves healing
- Elastic wrap in figure of 8 model improves ankle stability
- Hook & Loop closures ensure controlled compression and pain relief
- Double action design helps in prevention of inversion or eversion injuries



2508

Size: Small, Medium, Large, X-Large, XX-Large

SEGO[®] ANKLE SUPPORT

- Anatomically contoured, light weight and stretchable fabric
- Open heel design with single way stitch at the joint ensures free and easy movement
- For recurrent ligament and soft tissue injuries, Achilles tendon and pain due to arthritis



2501

Size: Small, Medium, Large, X-Large

SEGO[®] ANKLE BRACE

- Unique "Figure of 8" brace for firm ankle support
- Elastic wrap improves ankle stability and helps in the prevention of inversion or eversion injuries
- Ideal to treat soft tissue injuries, sprain and strain



2505

Size: 16",18",20",22",24"

SEGO® BREATH ANKLE BRACE

- Breathable elastic brace provides air circulation and prevents body odour
- Unique "Figure of 8" brace for firm ankle support
- Elastic wrap improves ankle stability and helps in the prevention of inversion or eversion injuries
- Ideal to treat soft tissue injuries, sprain and strain



1280

Size: Small, Medium, Large, X-Large

DYNA® ANKLE FOOT ORTHOSIS

- Made of polypropylene sheet which can be trimmed and moulded as per requirement
- Ideal for Flaccid Foot Drop and post operative stabilisation

Left & Right



1273

Size: Small, Medium, Large, X-Large

EASYCAST™ PNEUMATIC WALKER

- Light weight and durable shell provides optimal protection and comfort
- Non-skid rocker sole helps in pressure off-loading and promotes a natural gait
- Collateral air-cells offer pneumatic compression for total contact fit
- Inbuilt inflation and deflation pump allows custom inflation
- Toe bumper for protection

Available in 11" & 17" height



1275

Size: Small, Medium, Large, X-Large

DYNA® CAM WALKER

- Rocker Bottom promotes natural gait.
- Strong and rigid medio-lateral support improves immobilisation
- Light weight shell improves mobility
- Foam padding provides cushioning
- Hook and loop straps for easy application and removal

Available in 11" & 17" height





SURGICAL DRESSINGS AND BANDAGES



4001

Sizes:5,6,7,5,8,10 & 15cm X 4 & 4.5m stretched length

TOP'CREPE® COTTON CREPE BANDAGE B.P

- Made of cotton blended with rayon
- Special Z & S knitting gives 100% extensibility
- Highly porous for better air circulation
- Ideal for musculoskeletal injuries



4020

Sizes:8,10 & 15cm X 4m stretched length

TOP'GRIP® COTTON AND RUBBER ELASTIC BANDAGE B.P

- 300% extensibility: Offers graduated compression
- 100% elasticity: Reinforces the compression
- Conforms to body contours even after the swelling subsides
- Ideal for orthopaedic, phlebological and dermatologic problems



REF.

Sizes:A(4.5),B(6.25),C(6.75),D(7.5),E(8.75)
F(10)&G(12)cm X 1m &10m length

TUBIFIX® ELASTICATED TUBULAR BANDAGE B.P

- Tubular design, no tying and taping required
- Offers consistent compression
- Prevents fluid accumulation
- Ideal for difficult to dress joints
- Ideal for orthopaedic, plastic surgery and gynaecology conditions

4035(1m), 4037(10m)



REF.

Sizes:5,7,5,8,10,15cm X 4/6m stretched length

TOP'O PLAST® ELASTIC ADHESIVE BANDAGE B.P

- Long lasting adhesion
- Stretchable and has pores for air circulation
- For surgical, gynaecological, orthopaedic and general conditions



4027

4-LB™ MULTILAYER COMPRESSION BANDAGING SYSTEM (COMBI PACK)*

- Designed to deliver effective and sustained compression through out the wear time
- Fewer dressing changes and extended wear-time saves time and cost of treatment
- Contents: Plain Stockinette, Orthopaedic Bandage, Cotton Compression Bandage, Cotton and Rubber Elastic Bandage B.P., Cohesive Bandage



4110

Sizes:5,6,7.5,8,10,15,20cm width

PLAIN STOCKINETTE-UNBLEACHED

- Hygiene and convenient
- For fixation of wound dressings
- Available in 20 meters length



4125

Sizes:5cm,7.5cm,10cm& 15cm width X 3m length

TOP'BAN® ORTHOPAEDIC BANDAGE

- Excellent padding for plaster cases
- Conforms to body contours
- Resistant to plaster penetration
- Used prior to the application of Plaster of Paris Bandage



4100

Sizes: 5cm width X 1m length,7.5,10,15cm width X 1.5m length,5,7.5,10,15cm width X 10m length

TOP C NET® RIBBED COTTON STOCKINETTE B.P

- Made of cotton
- Ideal skin liner to be applied prior to prosthesis and plaster cases



4061

Sizes: Small, Medium, Large, X-Large

Thermafix[®] THERMAL SUPPORT

- Made of tubular knitted fabric for firm support
- Foam lining for support and comfort
- For post surgical/ injury support, arthritic conditions, sprain, strain and joint effusions

4061

Sizes: 1,2,3,4,5,6,7,8cm width X 25m stretched length

EASY NET[™] TUBULAR ELASTIC NET FIXATOR




PRODUCT FEATURES

- Allows quick application
- Hold wound dressings firmly in place
- Easy wound inspection
- Provides maximum aeration to wound
- Does not fray, all thread joints are knotted
- Available in different size, suitable for any part of the body
- Individual dispenser box for easy extraction of the bandage



ADVANCED DRESSINGS

SECUREMENT DRESSINGS



4165

Model: EC1 (Small, Medium, Large)

EASYFIX® Clear SEMI PERMEABLE ADHESIVE DRESSING

- Clear, breathable PU film provides a waterproof sterile barrier against microbes
- Flexible film conforms to body contours and stretches easily to prevent stress on the skin
- Skin- friendly adhesive is gentle to the patients skin and prevents accidental movement of catheters
- Extra strip provides additional securement
- Non -woven backing

Model No.	EC1		
Size	Small	Medium	Large
	4.1X6.5cm	5.5x6.5cm	6X8cm
Model Design			
Indication	IV Cannula Fixation		




4190

Size: Small, Medium, Large

EASYFIX® ELASTIC ADHESIVE BANDAGE B.P

- Convenient , saves nursing time
- Prevents infusion related phlebitis
- Stretches with body movements and allow air circulation
- Extra strip provides additional securement

Small	Medium	Large
3.5X6cm	5X6cm	6X8cm



4160

Size: EC35 (8.5cmX 10.5cm)

EASYFIX[®] Clear TRANSPARENT FILM DRESSING WITHOUT PAD (CATHETER DRESSING)

- Clear, breathable PU film provides a waterproof, sterile barrier against microbes
- Specially designed for challenging IV applications including the difficult to dress jugular sites
- Non Woven reinforcement with notched design improves catheter stability and dressing security
- Pre-cut sterile tape strips help no anchor hubs, lumens and infusion lines

Model No.	EC 35
Size	8.5cm X 10.5cm
Indications	Central line Catheterisation Dialysis Catheterisation



4160


Size: EC23 (6cmX 7cm), EC 24 (6cm X7cm)
EC26(10cm X 12cm)

EASYFIX[®] Clear SEMIPERMEABLE ADHESIVE DRESSING B.P

- Clear transparent PU film enables better and continuous monitoring of the insertion site
- Breathable material allows air circulation and reduces the risk of contamination
- The sterile film forms a protective barrier against bacteria virus and other external contaminants
- Water repelling property ensures that the site is dry
- Skin-friendly , hypoallergenic acrylic adhesive

Model No.	EC23	EC24	EC26
Size	6cm X7cm	6cmX7cm	10cmX12cm
Model Design			
Indications	IV Cannula Fixation	Laparoscopy Excisions Arthroscopy Carpal Tunnel Eye Cover	Hernioplasty Appendectomy Excisions Arthroscopy Epidural- Cathetherisation

PRE AND POST OPERATIVE DRESSINGS









4155

Size: EC 100(30cm X15cm), EC200(30cmX30cm),
EC 300(45cmX30cm), EC400(40cmX42cm),
EC500(55cmX45cm)

EASYFIX® Clear TRANSPARENT SURGICAL DRAPE

- Sterile film prevents bacterial regeneration and reduces the risk of wound contamination
- Adheres better to skin even during heavy retraction and irrigation
- Breathable, transparent, non-glare PU film conforms to the contour of the patient even on awkward operation sites

Type of drape	EC100	EC 200	EC 300	EC 400	EC 500
Size	30X15cm	30X30cm	45X30cm	40X42cm	55X45cm
Usage	Paediatric Surgery and Gynaecology	Paediatric and Hernia Surgery	Spinal and Abdominal Surgery	Hip, Knee replacement and Cardio Thoracic Surgery	Orthopaedic Thoracic and Vascular Surgery
Applications					







4153

Models: NW 70, NW71, NW72, NW73

STERIZONE® NW NON WOVEN ADHESIVE WOUND DRESSING

- Non-adherent wound contact layer prevents the dressing from sticking to the wound
- Highly absorbent pad absorbs exudate and minimises fluid strike through
- Non woven backing surface protect against extrinsic contamination, reducing the risk of infection
- Hypoallergenic acrylic adhesive

Sizes	NW70	NW71	NW72	NW73
Dressing Size	8.5X6cm	10cmX 15cm	20cmX10cm	25cmX10cm
Pad size	4cmX3cm	6cmX11cm	15cmX5cm	20cmX5cm
Model Design				

4170

EC 52, EC 54, EC 55, EC 56
EC 57, EC 58, EC 59, EC 60, EC 61, EC 63

EASYFIX[®] Clear POLYURETHANE WATERPROOF WOUND DRESSING



- Water - proof bacterial proof barrier
- Transparent PU film allows continuous monitoring of the wound area without removing the dressing
- Non- adherent absorbent pad absorbs excess wound fluid from the wound bed
- Doesnot stick to the wound and helps in faster healing

Sizes	EC 52	EC 56	EC 58	EC 61	EC 63
Dressing Size	5 X 7cm	9X10cm	15 X15 cm	9X25cm	9X35cm
Pad Size	2.5 X4cm	4.5X6.5cm	10X10cm	5X20cm	5X30cm
Usage	Laparoscopy Arthroscopy Peripheral Lines Skin Tear/ Abrasions Chronic Wounds	Appendix Removal Umbilical hernia Central lines Epidural Catheterisation Lacerations/ small trauma wounds Blisters/Small burns	Plastic Surgery: Graft Donor Site	OB/GYN C-Section Uterus Removal Breast Removal GENERAL SURGERY Gastrectomy Gallbladder Removal Hernia repair Kidney Removal Laminectomy	Open Heart Surgery Saphenous vein (Harvest site) Thoracotomy Total Hip Replacement Large Bowel Resection

Also Available

EC 54 (Dressing size: 5X 6.5cm, Pad Size: 2.5X 4cm)

EC 55 (Dressing size: 6X 10cm, Pad Size: 2.5X 6cm)

EC 57 (Dressing size: 8.5X 15.5cm, Pad Size: 4.5X 10cm)

EC 59 (Dressing size: 9X 15cm, Pad Size: 4.5X 10cm)

EC 60 (Dressing size: 9X 20cm, Pad Size: 4.5X 15cm)

4140


Size : 5X5cm, 10X 10cm

AgFix Foam[™] HYDROPHILIC FOAM SILVER PAD (NON-ADHESIVE)



ANTI ITCH
ANTI MICROBIAL
ANTI ODOUR

- Non- adherent wound contact layer ensure the dressing does not stick to the wound
- Anti- microbial silver layer kills 99.9% microbes and minimises the chances of infection
- Hydrophilic pad absorbs excess wound fluid reducing the risk of maceration
- Excellent pressure relief - ideal for diabetic , vascular and pressure ulcers



4150

Size: ST 80, ST 82, ST 83, ST86, ST 87, ST 88, ST 89, ST 90, ST 91, ST 93

STERIZONE™ SEMI PERMEABLE WATERPROOF PLASTIC WOUND DRESSING WITH SILVER

- Transparent film enables better and continuous monitoring of the area
- Breathable material provides a moist wound healing environment
- Absorbent pad absorbs excess wound fluid reducing the risk of maceration
- Silver layer kills 99.9% microbes and minimises the chances of infection
- Non- adherent wound contact layer prevents the dressing from sticking to the wound

Sizes	ST 80	ST 82	ST 86	ST 88	ST 91	ST 93
Dressing Size	45mm	5 X 7cm	9X10cm	15 X15 cm	9X25cm	9X35cm
Pad Size	24mm	2.5 X4cm	4.5X6.5cm	10X10cm	5X20cm	5X30cm
Usage	Laceration Abrasion minorcuts	Laparoscopy Arthroscopy Peripheral Lines Skin Tear/ Abrasions Chronic Wounds	Appendix Removal Umbilical hernia Central lines Epidural Catheterisation Lacerations/ small trauma wounds Blisters/Small burns	Plastic Surgery: Graft Donor Site	OB/GYN C-Section Uterus Removal Breast Removal GENERAL SURGERY Gastrectomy Gallbladder Removal Hernia repair Kidney Removal Laminectomy	Open Heart Surgery Saphenous vein (Harvest site) Thoracotomy Total Hip Replacement Large Bowel Resection

Also Available


- ST 83 (Dressing size: 5X 6.5cm, Pad Size: 2.5X 4cm)
- ST 84 (Dressing size: 6X 10cm, Pad Size: 2.5X 6cm)
- ST 87 (Dressing size: 8.5X 15.5cm, Pad Size: 4.5X 10cm)
- ST 89 (Dressing size: 9X 15cm, Pad Size: 4.5X 10cm)
- ST 90 (Dressing size: 9X 20cm, Pad Size: 4.5X 15cm)



COMPRESSION STOCKINGS

COMPREZON[®] MEDICAL COMPRESSION STOCKINGS

- Designed in line with the principles of graduated compression therapy
- Exert greatest pressure at the ankle, with compression decreasing gradually in the direction of the knee and thigh
- Stockings act as a layer of muscle by gently squeezing the vein walls, allowing the valves to close
- Restores blood flow to a normal state and aiding in overall circulation



Sizes: Small, Medium, Large, X-Large, XX-Large

COMPREZON[®] VARICOSE VEIN STOCKINGS CLASSIC


- Time tested and proven solution for venous disorders
- Offers medical compression of the highest level
- High tensile strength ensures durability and great elasticity

Styles: AD, AF, AG

Colours: Beige, White

Toe Styles: Peep hole, Open Toe, Closed Toe*

Class 1	AD	AF	AG	Class 2	AD	AF	AG
	2100	2101	2102		2110	2111	2112



2150


Sizes: Small, Medium, Large, X-Large, XX-Large

COMPREZON[®] VARICOSE VEIN STOCKINGS COTTON

- Made of 100% pure Egyptian cotton which is considered to be the world's finest
- Cotton absorbs moisture and remain gentle to the skin
- Ideal for all climates and ensures better patient compliance

Styles: AD, AF, AG

Class 1	AD	AF	AG	Class 2	AD	AF	AG
	2150	2151	2152		2160	2161	2162





Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

DVT-18[®] ANTI-EMBOLISM STOCKINGS

- Contour fitted graduated compression stockings
- Exert greatest pressure at the ankle, with compression decreasing gradually in the direction of the knee and thigh
- Stockings act as a layer of muscle by gently squeezing the vein walls, allowing the valves to close
- Restores blood flow to a normal state and aiding in overall circulation
- Prevents the danger of DVT

Styles	AD	AG
	2050	2052



2180

Sizes: Small, Medium, Large, X-Large, XX-Large

COMPREZON[®] THERAPEUTIC ENERGISING STOCKINGS

- Exerts 14-16 mmHg pressure at the ankle
- Over the calf stockings with graduated compression
- Gives massaging effect to the legs thereby energising the legs and reducing leg cramps, everyday leg swelling, mild oedema and tiredness



2180


Sizes: Small, Medium, Large, X-Large, XX-Large

easyfit[™] NURSE'S STOCKINGS

- Provides relief when standing for a very long time
- Graduated compression improves blood circulation
- Reduces mild leg swelling, leg pain and tiredness

STANDARD SIZE GUIDE(Lower Extremity)




Circumference in cm				
Sizes	Size identification colour	Circumference b	Circumference c	Circumference g
Small	Green	19-23	29-39	48-63
Medium	Blue	23-26	33-42	54-67
Large	Maroon	26-29	36-45	58-72
X-Large	Black	29-31	39-47	63-74
XX-Large	Red	31-34	44-52	74-86
XXX-Large	Light Blue	34-37	49-57	77-93



All length measurements have to be taken from point A

LENGTH MEASUREMENT IN CM

AD	AF	AG
37-42	55-70	70-80

Styles			Sizes						Toe Style	
			Small	Medium	Large	X-large	XX-Large	XXX-Large	Open Toe	
AD	AF	AG							Closed Toe	
									Peep Hole	

Size: Small, Medium, Large, X-Large, XX-Large



COMPREZON® LYMPHOEDEMA STOCKINGS/SLEEVES

- The elastic hosiery provides an intermittent massaging action thereby decongesting fluid build-up
- Aids in the upward flow of lymphatic fluid
- Ideal for post-mastectomy, post-traumatic or post-operative oedema and lymphoedema cases

Available Styles: AD, AF, AG, AGTR*/AGTL*, AT*, ATM*
(Lower extremity)
CG, CGH, AG, AGH (Sleeve)

* Available on request

AG	AGH	CG	CGH
2255	2265	2260	2270

STANDARD SIZE GUIDE (Lymphoedema Sleeve)


Available Styles




Armsleeve with hand



Armsleeve without hand




Armsleeve with hand + shoulder cap + belt



Armsleeve without hand + shoulder cap + belt

Available Sizes

Length C-G	standard 40-44cm, long 44-48cm				
	Sizes	Small	Medium	Large	X-Large
cA	18-21	20-23	22-25	24-27	26-29
cC	15-17	17-19	19-21	21-23	23-25
cE	23-29	26-32	29-35	32-38	35-41
cG	25-32	29-36	33-40	37-44	41-48





2195

Size: Small, Medium, Large, X-Large


VIBROX FLIGHT SOCKS

- Over the calf stockings with graduated compression
- Gently massages the calf muscles and improve venous circulation
- Energises legs and de-stress your feet
- Reduces leg swelling, pain and discomfort
- Ideal for flight and long journeys

Colours: Black, Blue



2219

Size: Small, Medium, Large, X-Large


VIBROX WORK SOCKS

- Over the calf stockings with graduated compression
- Provides relief when standing for a very long time
- Reduces mild leg swelling, leg pain and tiredness
- Energises legs and de-stress your feet
- Ideal for everyday wear

Colours: Black, Blue



TRACTION APPLIANCE



1504

One size fits most

DYNÁ[®] HOME CERVICAL TRACTION

- Components: Bucks apparatus, Chin and occipital pads with connecting straps, spreader bar and water bag
- Foam padded chin and occipital pads with over the door hanger
- For providing appropriate degree of traction to the cervical region



1502

One size fits most

DYNÁ[®] CERVICAL TRACTION

- Foam padded chin and occipital pads, connected with plastic cords to the spreader bar
- For providing appropriate level of traction at the cervical region both in sitting and lying position



1550

Size: Small, Medium, Large, X-Large, XX-Large

DYNÁ[®] PELVIC TRACTION

- For providing correct degree of traction at lumbo sacral region
- For soft tissue contractures and displaced vertebral fractures



1520

Size: Child, Adult

DYNÁ[®] LEG TRACTION

- Non- adhesive leg traction with non- slip foam lined stirrup attachment, straps and cords
- Used for both above and below knee after injury or surgery



1560

Size: Adult, Child

DYNA[®] LEG TRACTION APPLIANCE PREMIUM

- Consists of foam padding and cords for effective traction
- Offers even traction throughout the area of application
- Ideal for leg fracture, severe low back pain and when traction is required for a very long time



1540

Size: Small, Medium, Large, X-Large

DYNA[®] ANKLE TRACTION

- Made of high quality foam
- Can be ideally recommended for severe low back pain, hip/ sacral pain etc.



1570

Size: Child, Adult

SKIN TRAKT[®] ADHESIVE SKIN TRACTION KIT

- Consists of polyurethane foam lined stirrup attachment, Zinc Oxide Adhesive Bandage, straps and cords
- For severe low back pain and compound fractures



1508

Size: One size fits most

DYNA[®] BED TRACTION

- Used as an apparatus for suspending weights and fixing pulley belts to facilitate traction

1505

Size: One size fits most

DYNÁ TRACTION WEIGHT BAG



- To be filled with water for use in traction kit

A middle-aged man with a grey beard and glasses, wearing a light blue long-sleeved shirt, a black cap, and khaki trousers, is standing on a patterned carpet. He is using a four-wheeled walker with a black frame and copper-colored legs. The background shows a staircase with white railings and a building facade.

REHABILITATION AIDS



1942

Size: One size fits most

REHAID® ELBOW CRUTCHES WITH STRAP

- Made of light weight anodised Aluminium tube with telescopic height adjustment
- Padded elbow cuff provides comfort and total forearm support
- Additional hook & loop straps for stability
- PVC tip for better grip on the floor



1941

Size: One size fits most

REHAID® HINGED ELBOW CRUTCHES

- The hinged elbow cuff allows unrestricted movement of the wrist while providing firm support
- PVC tip for better grip on the floor
- Makes daily life activities easier



1938

Size: One size fits most

REHAID® CRUTCH STICK (SWAN NECK)

- Made of anodised Aluminium tube with telescopic height adjustment
- Swan neck shaped ergonomical handle
- PVC tip for better grip on the floor
- Adjustable height of 74cm to 94cm



1901

Sizes: Small, Medium, Large

REHAID® *Premium* ALUMINIUM AXILLARY CRUTCH

- With "C" shape Aluminium adjustable tube for better reinforcement
- Comfortable arm and palm padding made of EVA
- Heavy weight bearing capacity



1918

Size: One size fits most

REHAID[®] Premium T-SHAPE CRUTCH STICK

- Secured with ring locker for more safety
- T-shape design with premium look and feel



1932

Size: One size fits most

REHAID[®] CRUTCH STICK -TRIPOID

- Made of anodised Aluminium tube with telescopic height adjustment
- Three legged base for proper weight distribution and better support
- Swan neck handle with PVC tip for better grip on the floor
- Adjustable height of 78cm to 98cm



1943

Size: One size fits most

REHAID[®] CRUTCH STICK -QUADRIPOID

- Made of anodised Aluminium tube with telescopic height adjustment
- Four legged base for proper weight distribution and better support
- Swan neck handle with PVC tip for better grip on the floor
- Adjustable height of 78cm to 98cm



1910

Size: One size fits most

REHAID[®] Premium RECIPROCATING WALKING FRAME

- Can be used both as reciprocating and locked walking frame
- Easy to move forward without lifting the walker
- Ideal for gait training for people using prosthesis/orthosis, hemiplegic and geriatric patients
- Folds easily for convenient transporting



1912

Size: One size fits most

REHAID® WALKING FRAME WITH WHEELS

- Made of light weight anodised Aluminium tube with telescopic height adjustment
- Wheels help in making movement easy
- Ideal for gait training for people using prosthesis/orthosis, hemiplegic and geriatric patients
- Folds easily for convenient transporting
- Adjustable height of 77cm to 92cm



PODIATRY



1720

Size: Small, Medium, Large, X-Large, XX-Large
XXX-Large

SILICARE® SILICONE INSOLES(REMOVABLE)

- Biomechanically designed
- Made of skin friendly, medical grade silicone
- Soft, anti-shock, blue silicone zone incorporated at the heel and metatarsal region reduce skeletal shock and redistributes pressure
- Helps in maintaining the right arch
- Ideal for metatarsalgia, plantar fasciitis, tired and weak foot, joint, arch and back pain




1710

Size: Small, Medium, Large, X-Large

SILICARE® SILICONE HEEL CUSHION

- Biomechanically designed
- Made of skin friendly, medical grade silicone
- Gives protection and cushioning effect to the heel region
- Provides greater relief from external shock
- Ideal for heel and calcaneal pain, plantar fasciitis and bursitis of the heel

SILICARE® DIGITS AND FOREFOOT SUPPORT



METATARSAL CUSHION

For redistribution of pressure,painful metatarsalgia and collapsed transverse arch
Size :one size fits most



TOE SPREADER

For post Bunionectomy and Hallux Valgus
Sizes:Small,Medium,Large



BUNION SHIELD WITH RING

For Bunion
Size: One size fits most



OVAL DISC FOR CORNS

For corns and lesions between the toes
Size: One size fits most



METATARSAL PAD WITH ARMORED LOOP

For metatarsalgia, corns and calluses
Size: Medium, Large



STICK ON PADDING

For lesions, abrasions, cracking skin and heel blisters
Size: One size fits most




TOE SEPARATOR WITH SUPPORT

For overlapping toes, post forefoot surgery and interdigital lesions
Size: Small,Medium,Large



TOE CREST WITH LOOP

For hammer toes, claw toes and Morton's Neuroma
Size:Small,Medium,Large



OVAL DISC WITH HOLE FOR CORNS

For corns and painful lesions
Size: One size fits most



ROUND DISC WITH HOLE FOR CORNS

For interdigital soft corns
Size: One size fits most



DIGITAL RING

For corns and sores
Size: Small, Medium, Large



DIGITAL PROTECTOR

For soft corns and sores in the digits
Size: Small, Medium, Large



SELF ADHESIVE HEEL PROTECTOR

Ideal for shoe bites and irritation of achilles tendon
Size:Universal



TOE CAP

Nail injuries, digital ulcers and lesions at toe
Size: Small, Medium



3270

Sizes: Small, Medium, Large, X-Large

DYNA® ORTHOPAEDIC CAST FOOTWEAR

- Anti-skid rocker bottom sole supports forward movement
- For safe walking of patient with lower extremity cast and for cast protection
- For fixed claw toe, mid foot amputation or calcaneal ulcers



3280

Sizes: Small, Medium, Large, X-Large

DYNA® ORTHOPAEDIC SURGICAL FOOTWEAR

- Facilitate support to the healing foot after foot and ankle surgery
- Useful when foot cannot be fitted into normal footwear due to the presence of bandages and dressings



3290

Size: X-Small, Small, Medium, Large

DYNA® ANTERIOR ORTHO WEDGE FOOTWEAR

- Offloads pressure from forefoot to heel



3320

Size: X-Small, Small, Medium, Large

DYNA® POSTERIOR ORTHO WEDGE FOOTWEAR

- Offloads pressure from the heel to forefoot



MOTHER CARE



2700

Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

Cling[®] POST MATERNITY CORSET 25

- Made of premium quality elastic material
- 25cm wide corset covers the entire abdomen
- Supports weakened abdominal muscles post surgery/delivery



2720

Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

Cling[®] BREATH POST MATERNITY CORSET

- Made of premium quality breathable elastic material
- Special wide weave aids air circulation and prevents body odour
- 25cm wide corset covers the entire abdomen
- Supports weakened abdominal muscles post surgery/delivery



2730

Sizes: Small, Medium, Large, X-Large, XX-Large, XXX-Large

SEGO[®] ABDOMINAL CORSET PLAIN

- 20cm wide corset made of premium elastic material
- Fits the body contours and tightens abdominal muscles



2710

Sizes: Small, Medium, Large, X-Large, XX-Large
XXX-Large

SEGO® BREATH ABDOMINAL CORSET PLAIN

- 20cm wide corset made of breathable elastic material
- Fits the body contours and tightens abdominal muscles



2755

Size: Universal

SEGO® MATERNITY BINDER

- Provides excellent support to the lower abdomen during pregnancy
- Relieves added pressure on the spine caused due to weight gain
- Eases lower back pain and discomfort during pregnancy
- Can be adjusted to accommodate the abdomen during various stages of pregnancy



2734

Sizes: Universal

SEGO® ELASTIC ABDOMINAL CORSET (4 PANEL)

- Hook compatible elastic permits a closure at any point on the circumference of the binder
- Four panel design help in vertical size adjustments and reduce rolling
- Made of full elastic, soft material to provide comfort and flexibility
- Designed for use by both men and women



5101

Styles: Medi, Maxi

newmom™ DISPOSABLE MATERNITY PAD

- SAP (Super Absorbent Polymer) ensures maximum absorbency and less frequent pad changes
- Extra length for comfortable protection
- Stretches with body movements
- Prevents back and side flow ensuring maximum protection

Style	MEDIPAD	MAXIPAD
Absorbency Capacity	450ml	850ml



5155

Sizes: Small, Medium, Large, X-Large, XX-Large
XXX- Large

newmom™ DISPOSABLE PAD FIXATOR

- Breathable and non-constricting material- Reduces odour and ensures hygiene
- Soft and stretchable waist band- Prevents itching and pain caused by rubbing against the incisions
- Comfortably shaped for easy application and removal



REF.

Sizes: 60 X 90cm, 60 X 60cm

EASYPAD™ DISPOSABLE UNDERPAD

- Super Absorbent Polymer prevents re-wetting and protects skin from moisture
- Soft non-woven material prevents skin- irritation and bed- sores
- Sealed side prevents leakage
- Non- slip backing prevents moving of the pad

new mom™

∞ seamless



0820

Size: Small, Medium, Large, X-Large

new mom™ NURSING BRA

- Easy-to-open, double layered cups
- Wide, adjustable shoulder straps offer support and lift
- Wide back strap with 3 hook & eye closures for added support
- Adapts to the changing breast size at all times

Available : Black, Beige, White



0807

Size: Small, Medium, Large, X-Large

new mom™ C-SECTION PANTY

- Soft, non- constricting c - section panel doesn't impinge stitches
- Breathable panel aids in faster healing
- Light compression around the abdomen reduces swelling



0803

Size: Small, Medium, Large, X-Large

new mom™ MATERNITY SUPPORT PANTY

- Integrated support band below the tummy for reinforcement and lift
- Helps ease pressure on the lower back and belly
- Promotes proper posture



0801

Size: Small, Medium

new mom™ MATERNITY SUPPORT BELT

- Stretches to accommodate the growing tummy without being tight
- Helps ease pressure on the lower back & belly
- Promotes proper posture



0805

Size: Small, Medium

new mom™ MATERNITY HIPSTER

- Below-the-belly design comfortably supports the belly without compressing
- Stretchable elastic material progressively widens as the baby grows
- Breathable, soft, non-irritating material is ideal for delicate skin



0810

Size: Small, Medium

new mom™ COMFORT PANTY

- Seamless technology reduces chafing & irritation, provides all day comfort
- Stretchable fabric leaves 'no-marks' around the waist and thigh line
- Breathable material allows air circulation



SPORTS SUPPORTS



2550

Size: Small, Medium, Large, X-Large, XX-Large
XXX-Large

Olympian® KNEE SUPPORT

- Anatomically contoured cotton support provides firm and effective compression
- Ideal for various sports applications and day-to-day use
- For arthritis, chronic degenerative joint disease and irritation of the knee



2560

Size: Small, Medium, Large, X-Large, XX-Large

Olympian® KNEE SUPPORT LONG

- Anatomically contoured cotton support with increased length for extra support and effective compression
- Ideal for various sports applications and day-to-day use
- For arthritis, chronic degenerative joint disease and irritation of the knee



2510

Size: Small, Medium, Large, X-Large, XX-Large

Olympian® ANKLE SUPPORT

- Anatomically contoured cotton support provides firm and effective compression
- Ideal for various sports applications and day-to-day use
- For recurrent ligament & soft tissue injuries, Achilles tendon injuries, arthritis etc.



2901

Size: Small, Medium, Large, X-Large, XX-Large

Olympian® ELBOW SUPPORT

- Anatomically contoured cotton support provides firm and effective compression
- Ideal for various sports applications and day-to-day use
- For pain due to arthritis, weak & over stressed elbow and entrapment of median nerve in the elbow



2590

Size: Small, Medium, Large, X-Large, XX-Large

Olympian® THIGH SUPPORT

- Anatomically contoured cotton support provides firm and effective compression
- Ideal for various sports applications and day-to-day use
- For weak & over stressed quadriceps muscles, recurrent ligament & soft tissue injuries, sprain, strain and stiffness



2758

Size: Small, Medium, Large, X-Large, XX-Large

Olympian® ATHLETIC SUPPORT

- Ideal for protection and support during sports activities
- Designed to relieve fatigue and strain on the scrotal area



2924

Size: Small, Medium, Large, X-Large, XX-Large

Olympian® WRIST SUPPORT

- Unique "Pull and Wrap" design ensures proper fit
- Ideal for tendon injuries, sprain, strain, stiffness and arthritic conditions of the wrist
- Useful for various sports applications and day-to-day use



2932

Size: X-Small, Small, Medium, Large, X-Large, XX-Large, XXX-Large

Olympian® EXTENDED WRIST SUPPORT

- Knitted elastic support with extended length
- Ideal for tendon injuries, arthritis, sprain, strain and stiffness

2930

Size: X-Small, Small, Medium, Large, X-Large, XX-Large



Olympian® WRIST WRAP WITH THUMB SUPPORT

- Provides firm and effective support to the wrist joint and thumb
- For pain due to arthritis, sprain, strain and stiffness



0172

Size: Small, Medium, Large, X-Large, XX-Large

Lively™ ELBOW SUPPORT WITH STRAP

- Neoprene material offers soft tissue compression and therapeutic heat
- Retention which helps relieve elbow pain and promotes healing
- Contoured design provides needed support to counteract pain, oedema or inflammation
- Adjustable Hook & loop strap applies direct pressure on the tendons whilst the sleeve covers the whole elbow maximising the coverage



0170

Size: Universal

Lively™ ELBOW SUPPORT

- Soft neoprene support provides soothing warmth to muscles for improved blood circulation
- Supports, compresses and partially immobilizes the elbow joint following sprain or strain
- Contoured design provides needed support to counteract pain, oedema or inflammation



0178

Size: Universal

Lively™ TENNIS ELBOW BRACE SILICONE PAD

- Wraparound brace with adjustable fasteners for customized fit
- Silicone forearm pressure pad offers superior compression for faster relief
- Neoprene support provides warmth, reduces muscle fatigue and promotes blood circulation



0182

Size: Universal

Lively™ WRIST SUPPORT

- Soft neoprene support provides soothing warmth to the muscles for improved blood circulation
- Supports, compresses and partially immobilizes the wrist joint following sprain or strain
- Provides needed support to counteract pain, oedema or inflammation



0185

Size: Universal

Lively™ WRIST WRAP WITH THUMB SUPPORT

- Neoprene support provides warmth, reduces muscle fatigue and promotes blood circulation
- Offers therapeutic compression to weakened and strained wrist
- Prevent chances of injury to the wrist during sports or daily activities



0190

Size: Small, Medium, Large, X-Large

Lively™ SHOULDER SUPPORT

- Neoprene support provides warmth to the shoulder muscles and promotes blood circulation
- Provides optimal compression to the shoulder muscles and eases pain and swelling
- Strong grip on the shoulder with smooth and controlled assistance provides protection and support



0155

Size: Small, Medium, Large, X-Large, XX-Large

Lively™ THIGH SUPPORT

- Neoprene support provides warmth, reduces muscle fatigue and promotes blood circulation
- Provides optimal compression to the quadriceps and hamstring muscles and provide support during occupational and sporting activities
- Universal wraparound design enables easy application and customised adjustment



0158

Size: Small, Medium, Large, X-Large, XX-Large

Lively™ SHIN & CALF SUPPORT

- Neoprene support provides warmth, reduces muscle fatigue and promotes blood circulation
- Provides optimal compression to the shin and calf muscles and provide support during occupational and sporting activities
- Universal wraparound design enables easy application and customised adjustment



0146

Size: Small, Medium, Large, X-Large

Lively™ HINGED KNEE BRACE

- Neoprene support provides warmth, reduces muscle fatigue and promotes blood circulation
- The open patella design allows more flexibility around the knee cap and reduces pressure on the patella
- Soft circular cushion around the open patella prevents displacement of the knee cap and provides gentle massage
- Biaxial hinges allow natural flexion & free movement of the knee joint while eliminating unwanted gliding & twisting of the knee



0145

Size: Small, Medium, Large, X-Large

Lively™ KNEE SUPPORT OPEN PATELLA

- Neoprene material provides warming relief to the knee and aids in blood circulation
- The open patella design allows more flexibility around the knee cap and reduces pressure on the patella
- Soft circular cushion around the open patella prevents displacement of the knee cap and provides gentle massage



0140

Size: Universal

Lively™ KNEE SUPPORT

- Made of neoprene material which retains body heat and increases blood circulation
- Provides compression and support to the weak or injured knee
- Relieves pain due to strain or stiffness
- Open patella relieves pressure while soft buttress pad track the patella in position
- Wrap around design with opposing hook & loop straps for easy application and comfortable fit



0160

Size: Small, Medium, Large, X-Large

Lively™ ANKLE SUPPORT

- Soft neoprene support provides soothing warmth to the muscles for improved blood circulation
- Supports, compresses and partially immobilizes the ankle joint following sprain or strain
- Contoured design provides needed support to counteract pain, oedema or inflammation



0159

Size: Universal

Lively™ ANKLE WRAP

- Neoprene support provides warmth, reduces muscle fatigue and promotes blood circulation
- Innovative design contours to the shape of the affected area providing support to the injured muscles and joints for a speedy recovery
- Adjustable hook & loop straps provide targeted compression to reduce swelling
- Softer and stretchable material for better fit and adjustability




MISCELLANEOUS

REF.

Model: EC945, EC972

EASYTRODE™ ECG ELECTRODE

- Ag/AgCl electrode for excellent transmission
- Signal in 5 seconds after defibrillation and baseline in 2 seconds
- Impermeable foam backing with tab for easy peel off
- Skin friendly adhesive
- Recommended for resting ECG, TMT/Stress test, Holter ECG, In ICU, In operation theaters



5051

EASYGLIDE™ SURGICAL SKIN PREP BLADE

STERILE EO

- For Pre-Operative hair removal
- The product is sterile, individually packed for use in hospitals / clinics / ICU
- The design gives perfect angle for shaving hair
- The open design provides optimal view of shaving area





REF.

Style: Adult foot and ankle, Adult short leg, Adult hand
Adult short arm, Adult wide short arm, Adult long arm

DYNA[®] CAST PROTECTOR

- Reusable water proof plastic covering that keeps casts, bandages and dressings dry during bathing and showering
- Synthetic rubber diaphragm protects against water penetration while taking showers
- Protects casts, bandages, burns and abrasions
- Easy to apply and reusable



1701

Size: 6X6cm, 6X12cm, 3X15cm

SCAR TREAT[®] ADHESIVE GEL SHEET FOR SCARS

- Specially designed to soften, flatten, fade red dark or raised scars
- Works effectively on recent scars but also on older scars
- Bio-compatible and skin-friendly, will not cause allergy
- Easy to use, can be worn even on body contours and cut to fit the scar size



6040

easy pill® PILL ORGANISER

- An easy way for accurate medication
- Pill chambers for each day
- Individual slots for morning, noon, evening and night
- Pill chambers for the particular day can be easily detached and carried separately
- Braille letters to help visually challenged users



4146

Sizes: 5 X5cm, 10X10cm

PROZOLE™ ADHESIVE FELT PAD

- High resilient medical grade felt foam
- Offloads and redistributes pressure from the ulcer site
- Prevents abrasion to the wound margins
- Hypo allergenic skin friendly adhesive
- Allows direct fixation on skin

DISPLAY



Display Stand



Rev.No/Date:02B/19-03-2018



Dynamic Techno Medicals Pvt. Ltd

P.B. No.45 | Asokapuram | Aluva-683 101 | Kerala | India | Tel: +91 484 2837788

E-mail: export@dynamictechnomedicals.com | Website:www.dynamictechnomedicals.com